





2026 Patient Positioner Catalog

NinoMed.com



Safe-T-Secure® Family of Products





US Patents: 10512578, 11026854, 12023282, 12636213
European Patents: 2988717, 3545932, 3939557
Additional US and International Patents Pending

Our Original, All-In-One Trendelenburg patient positioner for robotic and laparoscopic surgery. It includes integrated handles to assist in moving the patient after intubation and an integrated perineal protective sheet to keep the OR table clean.

Safe -  - Secure®
by NinoMed, LLC



Integrated
Handles

Perineal Protective
Sheet

Better
GET A ^ GRIP!™

SIMPLE • EFFICIENT • EFFECTIVE

US Patents: 10512578, 11026854, 12023282, 12636213 | European Patents: 2988717, 3545932, 3939557 | Additional US and International Patents Pending



ORDER TODAY: (919) 818-3711 | NinoMed.com



Replace traditional solutions with the All-In-One Integrated Solution



VS

Traditional Solutions

- Egg crate mattress
- Foam for head
- Arm foam (2 per arm)
- Separate foam for hands
- Tape
- Gel pad
- Bean bag
- Arm boards
- Arm board disposable covers
- Shoulder blocks
- Linen or paper sheet
- Any additional straps
- 15-20 minutes for setup and removal

BENEFITS

UNIQUE WELL-DESIGNED CUSTOM SHAPE:

Designed by a laparoscopic and robotic surgeon to maximize efficiency and minimize setup time and effort.

SINGLE PATIENT USE WITH UNIQUE INTEGRATED INFECTION BARRIER:

Single use premium foam and impermeable perineal barrier designed to eliminate direct patient skin contact decreasing the risk of potential contamination.




EGG CRATE FOAM DESIGN:

Allows IV tubing and monitor wires to be tucked within the foam "valleys".

EASY AND EFFICIENT SETUP:

Fully integrated design includes all the foam, Velcro, and sheets necessary to secure a patient to the operating table. No need for arm boards.

Designed to completely stabilize patients in less than 1 minute, saving an average of 15-20 minutes of OR setup time.



Product Code

STSNMED1

- Vacuum sealed
- Efficient storage
- 5 per case

Safe-T-Secure® by NinoMed, LLC

Phone: 919-818-3711

Fax: 919-869-2012

Email: orders@ninomed.com



STS-BodyPad™ Positioning Kits



STS-BodyPad Regular™

Item #: STSBP1RG



Kit contains: 21" x 29" STS-Body Pad™, with integrated handles to assist with moving the pad plus patient down to the desired location, Lift Sheet, and Body Strap + Protective Foam.



STS-BodyPad Regular Plus™

Item #: STSBPRG2



Kit contains: 21" x 29" STS-Body Pad™, with integrated handles to assist with moving the pad plus patient down to the desired location, Lift Sheet, and Body Strap + Arm Protectors.



STS-BodyPad XL™

Item #: STSBP1XL



Kit contains: 21" x 40" STS-Body Pad™ - XL, with integrated handles to assist with moving the pad plus patient down to the desired location, Lift Sheet, and Body Strap + Protective Foam.



STS-BodyPad XL Plus™

Item #: STSBPXL2



Kit contains: 21" x 40" STS-Body Pad™ - XL, with integrated handles to assist with moving the pad plus patient down to the desired location, Lift Sheet, Body Strap + Protective Foam, plus Arm Protectors.

Pluma-Soft Head-to-Toe™ Pressure Redistribution Pad



Pluma-Soft™

Pressure Distribution Head-to-Toe™

Item #: STSH2T01



*Includes Padded Safety Strap
and Non-Woven Lift Sheet*

Take the guesswork out of pressure ulcer prevention and get covered head-to-toe. OR table pad designed for pressure ulcer prevention that spans the entire length of the OR Table. Includes non-woven lift sheet + body strap and protective foam.

Head and Neck Positioners

- Easy and consistent head and neck support during surgery
- Reliable pressure redistribution
- Option with slot to accommodate endotracheal tube
- Helps maintain proper alignment of the head, neck and spine during surgical procedures

Slotted Foam Head Positioner



Item #: SHPNM009,
SHPNM009-ECO

Description:
Slotted Foam Head Positioner
(8.5" x 8" x 4")

Order Qty: 24/case

Available in standard and
EcoEnhanced™ version

Multi-Ring Head Positioner



Item #: MRPNM012,
MRPNM012-ECO

Description:
Multi-Ring Head Positioner
(9" x 2")

Order Qty: 36/case

Available in standard and
EcoEnhanced™ version

Small Ring Head Positioner



Item #: SRDNM017,
SRDNM017-ECO

Description:
Small ring head positioner (7" x 2")

Order Qty: 72/case

Available in standard and
EcoEnhanced™ version

Adult Ring Head Positioner




Item #: FHDNM011,
FHDNM011-ECO

Description:
Ring Head Positioner (9" x 2")

Order Qty: 32/case

Available in standard and
EcoEnhanced™ version

Adult Head Positioner



Item #: AHPNM008,
AHPNM008-ECO

Description:
Adult Head Positioner (9" x 8" x 4.5")

Order Qty: 24/case

Available in standard and EcoEnhanced™ version

Cush-In It™ IV Wrap

Cush-In It™



CFNMED24
(4" x 1½")


Description:
Foam wrap designed to cushion a patient's IV tubing, monitor wires and accessories during surgery.

Order Qty: 6/pack; 50 packs/case

Upper Extremity Positioners

- Provides comfort and support for patient's arms and hands
- Designed to reduce the risk of pressure injury
- Flexible design for arms to be placed in various positions
- Ulnar nerve protector contains hook and loop strap

Arm Cradles



Item #: XACNM007,
XACNM007-ECO

Description:
Arm Cradles, Pair
(24" x 5" x 3")

Order Qty: 12 pairs/case

Available in standard and
EcoEnhanced™ version

Armboard Pad



Item #: APSNM006,
APSNM006-ECO

Description:
Armboard Pads, Pair
(20" x 8" x 2")

Order Qty: 12 pairs/case

Available in standard and
EcoEnhanced™ version

Ulnar Nerve Protector




Item #: UNPNM005,
UNPNM005-ECO

Description:
Ulnar Nerve Protectors, Pair

Order Qty: 36 pairs/case

Available in standard and
EcoEnhanced™ version

OR Table Pad



Item #: XRTPNM19,
XRTPNM19-ECO

Description:
OR table pad, convoluted
foam (2" x 20" x 72")

Order Qty: 12 pairs/case

Available in standard and
EcoEnhanced™ version

Lower Extremity Positioners

- Heel protectors are designed to reduce the risk of increased pressure during surgery
- Convoluted foam allows for air circulation and maximum comfort
- Features hook and loop strap

Heel Protector - Foam with Hook and Loop Strap



Item #: HPSNM013,
HPSNM013-ECO

Description:
Heel Protectors, Pair

Order Qty: 36 pairs/case

Available in standard and
EcoEnhanced™ version

Multi-Use Positioning Straps

- Single-use multi-purpose hook and loop straps
- Strong and dependable yet soft and comfortable
- Wide variety of applications: chest, body, and legs
- Stirrup straps available with and without slip ring

Positioning Body Straps - White



Item #: DBSNM001
Description:
Multi-Purpose White Positioning
Strap (60" x 4")
Order Qty: 12/case

Positioning Body Straps - Blue



Item #: XBSNM002
Description:
Positioning Strap, Blue
(72" x 5")
Order Qty: 20/case

Stirrup Strap with Slip Ring



Item #: SSWRNM03
Description:
Single-Use Stirrup Strap with Slip
Ring (19" x 3.5")
Order Qty: 240/case

Stirrup Strap without Slip Ring




Item #: XSNRNM04
Description:
Single-Use Stirrup Strap without
Slip Ring (19" x 3.5")
Order Qty: 240/case

Multi-Use Positioners

- Provides stability, comfort and elevation
- Can be used in a variety of procedures and patient applications
- Foam is supportive and comfortable for the patient

Large Foam Body Wedge



Custom sizes
available as well

Item #: LBWNM010,
LBWNM010-ECO
Description:
Large Body Wedge
(22" W x 11" D x 8" H)
Order Qty: 8/case
Available in standard and
EcoEnhanced™ version

Chest Roll



Item #: CRFNM016,
CRFNM016-ECO
Description:
Foam chest roll (16" x 5")
Order Qty: 12/case
Available in standard and
EcoEnhanced™ version



Safe-T-Secure® Family of Products.

Simple, effective and efficient solutions
for all your positioning needs.

Our mission is to provide innovative and useful medical products at affordable prices. NinoMed was founded by an academic surgeon with a passion for patient safety and operating room efficiency. We start by first understanding the challenges hospitals face, and then create products and solutions that improve patient care. NinoMed's positioners are designed to improve OR setup, efficiency, patient comfort and safety while reducing costs.

Patient Positioners Catalog
Additional and custom products
available on request.

Copyright © 2026 by NinoMed, LLC.
All rights reserved.

ORDER TODAY: (919) 818-3711 | NinoMed.com

US Patents: 10512578, 11026854, 12023282, 12636213 | European Patents: 2988717, 3545932, 3939557 | Additional US and International Patents Pending



NinoMed-Catalog-0005 7/10/26