





2023 Patient Positioner Catalog

[NinoMed.com](http://NinoMed.com)



# Smart Patient Solutions That Work





US Patent: 10512578, European Patent: 2988717  
Additional Patents Pending

THE ORIGINAL ALL-IN-ONE TRENDLENBURG PATIENT POSITIONING SOLUTION

Safe -  - Secure<sup>®</sup>  
by NinoMed, LLC

*A complete Trendelenburg  
positioning solution for all  
laparoscopic & robotic surgery*



*Better*  
**GET A ^ GRIP!™**

SIMPLE • EFFICIENT • EFFECTIVE

US Patent: 10512578  
European Patent: 2988717  
Additional Patents Pending



  
**NINO MED<sup>®</sup>**

## Replace Traditional Solutions with *The Original All-In-One Integrated Solution*



VS

### Traditional Solutions

- Egg crate mattress
- Foam for head
- Arm foam (2 per arm)
- Separate foam for hands
- Tape
- Gel pad
- Bean bag
- Arm boards
- Arm board disposable covers
- Shoulder blocks
- Linen or paper sheet
- Any additional straps
- 15-20 minutes for setup and removal

### Benefits:

#### UNIQUE WELL-DESIGNED CUSTOM SHAPE:

Designed by a laparoscopic and robotic surgeon to maximize efficiency and minimize setup time and effort.

#### SINGLE PATIENT USE WITH UNIQUE INTEGRATED INFECTION BARRIER:

Single use premium foam and impermeable perineal barrier designed to eliminate direct patient skin contact decreasing the risk of potential contamination.



#### EGG CRATE FOAM DESIGN:

Allows IV tubing and monitor wires to be tucked within the foam "valleys".

#### EASY AND EFFICIENT SETUP:

Fully integrated design includes all the foam, Velcro, and sheets necessary to secure a patient to the operating table. No need for arm boards.

Designed to completely stabilize patients in less than 1 minute, saving an average of 15-20 minutes of OR setup time.



#### Product Code

STSNMED1

- Vacuum sealed
- Efficient storage
- 5 per case





www.ninomed.com

US Patent: 10512578 European Patent: 2988717 Additional Patents Pending

EVERYTHING YOU LOVE ABOUT SAFE-T-SECURE® PLUS A FLEXIBLE BODY STRAP

# SIS-FLEX™

by NinoMed, LLC

*A flexible and complete positioning solution for both cases in Trendelenburg and cardiothoracic surgery*



Provides optimal access for cardiothoracic surgery

Uniquely designed for:

- Maximum flexibility
- Efficient and safe arm tucking
- Use in a wide variety of surgeries
- Pressure relief
- Cardiothoracic and abdominal surgery with its flexible body strap

SIMPLE • EFFICIENT • EFFECTIVE ... AND FLEXIBLE

US Patent 10512578  
European Patent: 2988717  
Additional Patents Pending



**NINO/MED®**

## STS-Flex™ is designed to:

- Easily position a patient on any surgical table.
- Secure a patient in a wide variety of surgeries with its flexible and adjustable safety strap.
- Reduce the risk of pressure ulcer formation by maximizing pressure redistribution.
- Reduce hospital acquired infections with its integrated perineal infectious barrier sheet.
- Keep the OR table clean and reduce turnover time with protective perineal barrier.
- Protect a patient's arms and hands with integrated foam wings that extend below the OR table.



## Benefits:

**UNIQUE WELL-DESIGNED CUSTOM SHAPE:** Designed by a laparoscopic and robotic surgeon to maximize efficiency and minimize setup time and effort.

**FLEXIBLE BODY STRAP:** Everything you love about Safe-T-Secure® plus a flexible body strap. The adjustable padded body strap allows for STS-Flex™ to be used in a wide variety of surgeries including cardiothoracic, gynecology, urology, colorectal and general surgery.

**SINGLE PATIENT USE WITH UNIQUE INTEGRATED INFECTION BARRIER:** Single use proprietary foam and impermeable perineal barrier designed to eliminate direct patient skin contact with OR table decreasing the risk of potential contamination.

**EGG CRATE FOAM DESIGN:** Allows IV tubing and monitor wires to be tucked within the foam "valleys".

**EASY AND EFFICIENT SETUP:** Integrated design includes all the foam, Velcro, and sheets necessary to secure a patient to the operating table. No need for arm boards.



### Product Code

STS-FLEX1

- Vacuum sealed
- Efficient storage
- 5 per case

## STS-Head-to-Toe™



Includes Padded Safety Strap  
and Non-Woven Lift Sheet



### STS-Head-To-Toe™

Item #: STSH2TO1

**Description:**

Take the guesswork out of pressure ulcer prevention and get covered head-to-toe. OR table pad designed for pressure ulcer prevention. Includes non-woven lift sheet + body strap/foam. (20" x 76" x 3.6")

Order Qty: 5/case

## STS-Flex-Extend™



### STS-Flex-Extend™

Item #: STSFLEXE

**Description:**

Full length OR table pad with integrated arm tucking. It has a flexible and adjustable padded body strap that can be attached anywhere on the OR table.

Order Qty: 5/case

## STS-BodyPad™



Includes Padded Safety Strap  
and Non-Woven Lift Sheet



### STS-BodyPad™

Item #: STSHAB01

**Description:**

OR table pad designed for lithotomy. Has perineal protective sheet. Includes non-woven lift sheet + body strap/foam. (20" x 40" x 3.6")

Order Qty: 5/case



## Head and Neck Positioners

- Easy and consistent head and neck support during surgery
- Reliable pressure redistribution
- Option with slot to accommodate endotracheal tube
- Helps maintain proper alignment of the head, neck and spine during surgical procedures

### Slotted Foam Head Positioner



Item #: SHPNM009,  
SHPNM009-ECO

Description:  
Slotted Foam Head Positioner  
(8.5" x 8" x 4")

Order Qty: 24/case

Available in standard and  
EcoEnhanced™ version

### Multi-Ring Head Positioner



Item #: MRPNM012,  
MRPNM012-ECO

Description:  
Multi-Ring Head Positioner  
(9" x 2")

Order Qty: 36/case

Available in standard and  
EcoEnhanced™ version

### Small Ring Head Positioner



Item #: SRDNM017,  
SRDNM017-ECO

Description:  
Small ring head positioner (7" x 2")

Order Qty: 72/case

Available in standard and  
EcoEnhanced™ version

### Adult Ring Head Positioner




Item #: FHDNM011,  
FHDNM011-ECO

Description:  
Ring Head Positioner (9" x 2")

Order Qty: 32/case

Available in standard and  
EcoEnhanced™ version

### Adult Head Positioner



Item #: AHPNM008,  
AHPNM008-ECO

Description:  
Adult Head Positioner (9" x 8" x 4.5")

Order Qty: 24/case

Available in standard and EcoEnhanced™ version

## Cush-In It™ IV Wrap




Item #s: CFNMEDV1  
(variety pack, 2 of each size)


### Cush-In It™



Small  
CFNMED22  
(2" x 1 1/4")



Medium  
CFNMED23  
(3" x 1 3/8")



Large  
CFNMED24  
(4" x 1 1/2")

#### Description:


Foam wrap designed to cushion a patient's IV tubing, monitor wires and accessories during surgery. Available in (S) 2" x 1 1/4", (M) 3" x 1 3/8", and (L) 4" x 1 1/2".

Order Qty: 6/pack; 50 packs/case

## Upper Extremity Positioners

- Provides comfort and support for patient's arms and hands
- Designed to reduce the risk of pressure injury
- Flexible design for arms to be placed in various positions
- Ulnar nerve protector contains hook and loop strap

### Arm Cradles



Item #: XACNM007,  
XACNM007-ECO

Description:  
Arm Cradles, Pair  
(24" x 5" x 3")

Order Qty: 12 pairs/case

Available in standard and  
EcoEnhanced™ version

### Armboard Pad



Item #: APSNM006,  
APSNM006-ECO

Description:  
Armboard Pads, Pair  
(20" x 8" x 2")

Order Qty: 12 pairs/case

Available in standard and  
EcoEnhanced™ version

### Ulnar Nerve Protector




Item #: UNPNM005,  
UNPNM005-ECO

Description:  
Ulnar Nerve Protectors, Pair

Order Qty: 36 pairs/case

Available in standard and  
EcoEnhanced™ version

### OR Table Pad



Item #: XRTPNM19,  
XRTPNM19-ECO

Description:  
OR table pad, convoluted  
foam (2" x 20" x 72")

Order Qty: 12 pairs/case

Available in standard and  
EcoEnhanced™ version

## Lower Extremity Positioners

- Heel protectors are designed to reduce the risk of increased pressure during surgery
- Convoluted foam allows for air circulation and maximum comfort
- Features hook and loop strap

### Heel Protector - Foam with Hook and Loop Strap



Item #: HPSNM013,  
HPSNM013-ECO

Description:  
Heel Protectors, Pair

Order Qty: 36 pairs/case

Available in standard and  
EcoEnhanced™ version



## Multi-Use Positioning Straps

- Single-use multi-purpose hook and loop straps
- Strong and dependable yet soft and comfortable
- Wide variety of applications: chest, body, and legs
- Stirrup straps available with and without slip ring

### Positioning Body Straps - White



Item #: DBSNM001  
Description:  
Multi-Purpose White Positioning  
Strap (60" x 4")  
Order Qty: 12/case

### Positioning Body Straps - Blue



Item #: XBSNM002  
Description:  
Positioning Strap, Blue  
(72" x 5")  
Order Qty: 20/case

### Stirrup Strap with Slip Ring



Item #: SSWRNM03  
Description:  
Single-Use Stirrup Strap with Slip  
Ring (19" x 3.5")  
Order Qty: 240/case

### Stirrup Strap without Slip Ring




Item #: XSNRNM04  
Description:  
Single-Use Stirrup Strap without  
Slip Ring (19" x 3.5")  
Order Qty: 240/case

## Multi-Use Positioners

- Provides stability, comfort and elevation
- Can be used in a variety of procedures and patient applications
- Foam is supportive and comfortable for the patient

### Large Foam Body Wedge



Custom sizes  
available as well

Item #: LBWNM010,  
LBWNM010-ECO  
Description:  
Large Body Wedge  
(22" W x 11" D x 8" H)  
Order Qty: 8/case  
Available in standard and  
EcoEnhanced™ version

### Chest Roll



Item #: CRFNM016,  
CRFNM016-ECO  
Description:  
Foam chest roll (16" x 5")  
Order Qty: 12/case  
Available in standard and  
EcoEnhanced™ version

## 2 OPTIONS, 1 GREAT SOLUTION

### Safe-T-Secure®



- The original, All-In-One Trendelenburg patient positioning solution
- Efficient and safe arm tucking
- Designed to improve OR efficiency and to safely and effectively position patients
- Integrated perineal infectious barrier designed to reduce hospital acquired infections



### Safe-T-Secure®

Item #: STSNMED1

Description:

Original, All-In-One Trendelenburg patient positioning solution for robotic and laparoscopic surgery

Order Qty: 5/case



US Patent: 10512578, European Patent: 2988717  
Additional Patents Pending

### STS-Flex™



- Everything you love about Safe-T-Secure® plus a flexible body strap
- Easily position a patient on any surgical table
- Features adjustable padded body strap
- Flexible body strap for both cardiothoracic and abdominal surgery



### STS-Flex™

Item #: STSFLEX1

Description:

A flexible and complete positioning solution for both cases in Trendelenburg and cardiothoracic surgery

Order Qty: 5/case



By NinoMed, LLC  
www.ninomed.com  
US Patent: 10512578, European Patent: 2988717, Additional Patents Pending

## INDEX

Safe-T-Secure® The Original All-In-One Trendelenburg Positioner .....	2
STS-Flex™ Trendelenburg and Cardiothoracic Positioner .....	4
ORDERING INFORMATION	
STS Head-To-Toe™ .....	6
STS-Flex-Extend™ .....	6
STS-BodyPad™ .....	6
Head and Neck Positioners .....	7
Cush-In It™ IV Wrap .....	7
Upper Extremity Positioners .....	8
Lower Extremity Positioners .....	8
Multi-Use Positioning Straps .....	9
Multi-Use Positioners .....	9
Safe-T-Secure® The Original All-In-One Trendelenburg Positioner .....	10
STS-Flex™ Trendelenburg and Cardiothoracic Positioner .....	10



### Simple, Effective and Efficient Solutions For All Your Positioning Needs

Patient Positioners Catalog  
Additional and custom products available on request.

Copyright © 2023 by NinoMed, LLC. All rights reserved.



# NINO/MED®

Our mission is to provide innovative and useful medical products at affordable prices. NinoMed was founded by an academic surgeon with a passion for patient safety and operating room efficiency. We start by first understanding the challenges hospitals face, and then create products and solutions that improve patient care. NinoMed's positioners are designed to improve OR setup, efficiency, patient comfort and safety while reducing costs.

**ORDER TODAY: (919) 818-3711 | [NinoMed.com](http://NinoMed.com)**

